

DV8 Steerable Aluminum Knee Walker



HCPCS: E0118









FIGURE A

FIGURE B

FIGURE C

FIGURE D







FIGURE F

FIGURE G

- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- Deluxe dual-braking system
- Sturdy dual frame adds stability and durability
- New basket allows the user to transport personal items (Figure D)
- Contoured leg pad provides extra balance
- New leg pad cover provides added comfort (Figure E)
- New tiller lock provides added safety and security
- Tiller can be easily folded for storage and transportation
- Foldable steering column allows for easy storage
- Tool-free height adjustment of seat and handle (Figure A B)
- 8" casters are ideal for indoor/outdoor use
- Silver vein finish (Figure C) is durable and easy-to-maintain

PRODUCTS	Description
790	Steerable Knee Walker,1/ea
SPECIFICATIONS	Description
Carton Shipping Weight	26 lbs.
Handle Height	33" - 37.5"
Knee Pad Height	17.5" - 21.5"
Unit Dimensions	17" (W) x 32" (D)
Warranty	Limited Lifetime, Lifetime Warranty on Brake Cable
Weight	19 lbs.
Weight Capacity	300 lbs
PARTS	Description
9503H79015-N	KneePadHeightKnob790New,1ea, CS
9503H79016	Left Wheel Lock,790,1ea, CS
9503H79017	Right Wheel Lock,790,1ea, CS
9503H79018	Caster Complete Assy,790,1ea, CS
9503H79019	Caster Axle,790,1ea, CS
9503H79020	Axle Set Screws,790,4/pkg, CS
9503H79021	Wheel Bearings,790,PhoenixHD3,1/pr, CS
9503H79001	Handle Bar,790,1ea, CS
9503H79002	Hand Grips,790,1Pair, CS



PARTS	Description
9503H79003	Left Brake Complete,790,1ea, CS
9503H79004	Right Brake Complete,790,1ea, CS
9503H79005	Handle Bar Column,790,1ea, CS
9503H79006	Handle Bar Height Adj Lock,790,1ea, CS
9503H79007	Steering Column,790,1ea, CS
9503H79008	Brake Cable Guide,790,1ea, CS
9503H79011	Steering Fork,790,1ea, CS
9503H79012	Column To Fork Locking Bolt,790,1ea, CS
9503H79013	Fork Bearings,790,1pair, CS
9503H79014	Knee Pad,790,1ea, CS
9503H79023	Brake Cable & Hardware Only for 790, 1ea, CS