DIMENSIONS	FDA, TÜV, CE Approved	
Trolley mode:	H 28"	71 cm
	L 15.3"	38.9 cm
	W 16.5"	41.9 cm
Drive mode:	H 35"	89 cm
	L 47"	120 cm
	W 22"	56 cm

WEIGHT

Front half:	26.4 lbs	12 kg
Rear half:	37.5 lbs	17 kg

29 kg

		(W/o patte
Battery:	4.1 lbs	1.85 kg

64 lbs

Battery: 4.1 lbs 1.85 kg

WEIGHT CAPACITY

Assembled:

Maximum	220 lbs	100 kg
user weight:	220 103	

BATTERY

48 Volts Lithium-ion, more than 10 miles of continuous driving, only 4 hours on average for full charge, charge on-board and off-board

SPEED

Most countries, 4 mph/ 6.4 km per hour*
* According to the local regulation



Separates into 2 parts in seconds



Easy to wheel as a trolley



Compact to full size in 10 sec



The best folding scooter in the world