

## Fly-Lite Aluminum Transport Chair

HCPCS: E1038



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Weighs only 16.2 lbs without footrests 18.9 lbs with footrests
- Aluminum frame is lightweight and strong
- Comes with carry pocket on backrest
- Comes with seat belt for added safety
- Attractive plaid style nylon upholstery (Figure A)
- Comes with handy carry pocket on back (Figure B)
- Tool-free, height-adjustable swing-away footrests come standard
- Optional leg rests available (Model# STDELR-TF)
- 8" rear wheel 6" casters in front
- Back folds down for storage and transport and features deluxe back release (Figure C)
- Comes with aluminum rear wheel locks
- Convenient cup holder can easily attach to chair Sold separately (Model# STDS1040S) (Figure D)
- Weight Capacity: 300 lbs

### PRODUCTS

#### Description

DFL19-BL	Fly-LiteTransportChair,19"Bl,1/cs
DFL19-BLK	Fly-LiteTransportChair,19"Blk,1/cs
DFL19-RD	Fly-LiteTransportChair,19 Red,1/cs

### SPECIFICATIONS

#### Description

Back Height	18"
Carton Shipping Weight	23.8 lbs
Overall Height	36"
Overall Length w/ Riggings	33"
Overall Width	22" (open)   9" (closed)
Seat Dimensions	19" (W) x 15.5" (D)
Seat To Floor Height	18"
Weight Capacity	300 lbs
Weight With Front Riggings	18.9 lbs
Weight Without Front Riggings	16.2 lbs

### ACCESSORIES

#### Description

STDS1040S	Cup Holder, Swivel,1BX,6/CS, CS
-----------	---------------------------------

### PARTS

#### Description

STDS2A1807	SEAT RAIL PAD,DFL,NITRO,1/EA, CS
STDS2A1803	Cross Frame Assist (Blue) DFL,1/EA, CS

PARTS	Description
STDS2A184R	BRAKE ASSY RIGHT,DFL,1/EA, CS
STDS2A184L	BRAKE ASSY LEFT,DFL,1/EA, CS
STDS2A1810	FRONT WHEEL,DFL,TS19,1/ea, CS
STDS2A1806	FRONT AXLE BOLT,DFL,1 PR, CS
STDS2A1811	REAR AXLE BOLT,DFL,1 PR, CS
STDS2A1801	REAR HUB CAP,DFL,1 PR, CS
STDS2A1800	Rear Wheel, DFL, TS19,1/ea, CS
TRB-1	BackSupportReleaseLever./Trans,1/ea, CS
STDS2A181L	HANDLE BAR,LEFT,BLACK,DFL,1/EA, CS
STDS2A182L	HANDLE BAR,LEFT,BLUE,DFL,1/EA, CS
STDS2A1891	Back Upholstery, Black/DFL,1/ea, CS
STDS2A1892	Back Upholstery, SilverPlaid,DFL,1/EA, CS
STDS2A1893	BACK UPHOLSTERY Blue,DFL,1/EA, CS
STDS001014	Folding Cane Hindge DFL,1/ea, CS
STDS2A1824	Arm Pad,Black,1/ea, CS
STDS2A1804	Arm Pad,Black Plaid,DFL,1ea, CS
STDS2A1814	Arm Pad,Blue Plaid,DFL,1ea, CS
DFLSFBK	DFL Footrests, Black, 1pr/cs, CS