

The Jazzy Select 6® 2.0 offers a truly unique ride experience with a sporty design and superior performance. Pride's patented Active-Trac® Suspension combines with Mid-Wheel 6® drive design for a smooth, comfortable and stable ride, both indoors and out. The depth-adjustable high-back seat folds flat for greater convenience, while offering maximum comfort and support.

THIS FOOT PLATFORM OFFERS SUPERIOR ADJUSTABILITY FOR AMPLE FOOT ROOM.









THIS CONVENIENT FOLDING SEAT IS DEPTH ADJUSTABLE FOR MAXIMUM COMFORT AND OFFERS REPLACEABLE BACK AND SEAT COVERS.

AMERICA'S FAVORITE POWER CHAIR

Learn more online.









gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed ¹	Up to 4.3 mph
Ground Clearance ²	2.5" at footplate, 1.77" at motor bracket
Turning Radius ²	23.76"
Overall Length ²	43.27"
Overall Width	23.27"
Drive Wheels	10" flat free
Front Wheels	6" flat free
Rear Wheels	6" flat free
Suspension	Active-Trac® (ATX) with Mid-Wheel 6® Technology
Drivetrain	Two-motor, Mid-Wheel 6® drive
Braking System	"Intelligent Braking" electronic regenerative, disc park brake
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3.5 amp
Per-Charge Range ^{2,4}	Up to 19 miles
Battery Requirements ^{3,5}	(2) 12 volt, deep cycle, U-1
Battery Weight	24.5 lbs. (U-1)
Base Weight	170 lbs.without batteries
Standard Seat Weight	50 lbs. (comfort, high-back)
Weight of Heaviest Piece	114 lbs.
Max. Seat Dimensions	20" x 20" (high-back)
Seat-to-Floor Range	21.5" x 23.25"
Warranty	5-year limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Features

- Active-Trac® Suspension combines with Mid-Wheel 6® Technology for a smooth, stable ride over uneven terrain
- Convenient, limited recline, high-back black vinyl folding seat is depth adjustable for maximum comfort and includes replaceable back and seat covers available in sizes: 16" x 16"-18", 18" x 18"-20" and 20" x 18"-20
- · Height-adjustable, flip-up foot platform offers ample foot room for maximum comfort
- · Black, flat-free, non-marking tires saves your floors
- · In-line joystick mount is equipped with performancematched electronics for easy operation
- Equipped with 60" lap belt for security and peace of mind

Optional Seating*

- Limited recline, high-back black vinyl solid seat pan with depth adjustment available in sizes: 16" x 16"-18" or 18" x 18"-20"
- · Deluxe contour, limited recline, high back with solid seat pan, available sizes: 20" x 18" or 20" x 20"
- Power elevating seat (VR2 6-key required)

Options

- Elevating leg rests
- Swing-away footrests
- Heel loops
- · Swing-away joystick
- Residual limb support

Accessories

- Cane/crutch holder
- Cup holder
- · Flag holder
- Lap belt (50", 70", 80", 90")

can be subject to a variance of (+ or - 3%).

5 Battery weight may vary based on manufacturer.

3 AGM or gel cell type required.

- Military patches
- · Oxygen tank holder
- · Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

Model

Jazzy Select 6 2.0 | Code: K0822/K0823

system performance. Test conducted at maximum weight capacity.

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/-2 Due to manufacturing tolerances and continual product improvement, this specification

Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards.

Results derived from theoretical calculations based on battery specifications and drive

Colors:



Blue



Ruby Red





- Sapphire
- - Pearl White
- - Black Onyx

*Additional charges may apply.











