drive DeVilbiss

ZooMe Flex Folding Travel Scooter



HCPCS:





FIGURE C



FIGURE D

FIGURE A FIGURE B



- Ultra-lightweight for superior portability
- Quick and Easy "Flex" folding design for transport and storage in one easy motion
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests and fold down back
- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller angle places the scooter's controls in the ideal position
- Anti-tip wheels for additional safety

PRODUCTS	Description
FLEX	ZooMe Flex Folding Travel Scooter, Blue
SPECIFICATIONS	Description
Anti-Tip Wheels	2" x .5"
Base Weight	53.6 lbs.
Batteries (pair)	25.9V x 10AH Lithium
Battery Charger	2A Offboard
Battery Weight (pair)	5 lbs.
Brakes	Electromagnetic
Climbing Angle	6°
Controller	PG S Drive 45A
Estimated Range	13 miles
Freewheel Mode	No
Front Wheels	6" x 1.5"
Ground Clearance	2"
HC-FEET	4" (L) x 3" (W) x 0.5" (H)
HC-HAND	5" (L) x 3" (W) x 1.25" (H)
Heaviest Peice Weight	53.6 lbs
Motor	24V x 120W
Overall Dimensions	37" (L) x 18" (W) x 36" (H)
Rear Wheels	7" x 2.35"
Seat Dimensions	16.75" (W) x 15" (D)
Top Speed	4 mph
Total Weight (excl. hand control)	58.6 lbs.
Turning Radius	53"

SPECIFICATIONS	Description
Warranty	24 months on Electronics, 12 months on Batteries, Limited Lifetime on Frame
Weight Capacity	275 lbs
Weight Capacity (13046 13046-HD)	275 lbs.
PARTS	Description
FLEX-01	Upper Panel Assembly,Blue,ZooMeFlex,1/ea, CS
FLEX-02	Throttle Assembly, ZooMe Flex, 1/ea, CS
FLEX-03	Gas Spring Assembly, ZooMe Flex, 1/ea, CS
FLEX-04	FrontReflectorAssmbly,Blue,ZooMeFlex,1ea, CS
FLEX-05	Throttle Potentiometer, ZooMe Flex, 1/e, CS
FLEX-06	Quick Release Clamp, ZooMe Flex, 1/ea, CS
FLEX-07	Front Wheel Assembly, ZooMe Flex, 1/ea, CS
FLEX-08	Steering Knuckle, Left, ZooMe Flex, 1/e, CS
FLEX-09	Steering Knuckle,Right, ZooMe Flex, 1/e, CS
FLEX-10	Hub Cap, ZooMe Flex, 1/ea, CS
FLEX-11	Wheel / Stem Bearings, ZooMe Flex, 1/pr, CS
FLEX-12	Front Cover, Blue, ZooMe Flex, 1/ea, CS
FLEX-13	Tie Rod W/Ends, ZooMe Flex, 1/ea, CS
FLEX-14	Handle Connector, ZooMe Flex, 1/ea, CS
FLEX-15	Steering Connector, ZooMe Flex, 1/ea, CS
FLEX-16	BatteryCase W/O Batteries,ZooMe Flex,1e, CS
FLEX-17	Front Pedal Cover, Blue,ZooMe Flex, 1/ea, CS
FLEX-18	RotatingGearAssembly,Front,ZooMeFlex,1e, CS
FLEX-20	ReleaseCableW/Hrdwr&Lever,ZooMeFlex,1/ea, CS
FLEX-19	Support Rod, Front, ZooMe Flex, 1/ea, CS