

Blue Streak Wheelchair

HCPCS: K0001-1,K0001



FIGURE A



FIGURE B

- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- Solid rubber tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap
- Comes with swing-away footrests or elevating leg rests
- Weight Capacity: 250 lbs

PRODUCTS

Description

BLS16FBD-ELR	BlueStreak16"D/L FlipBackArms,Elev,1/ea
BLS16FBD-SF	BlueStreak16"D/L FlipBackArms,SwgFt,1/ea
BLS18FBD-ELR	BlueStreak18"D/L FlipBackArms,Elev,1/cs
BLS18FBD-SF	BlueStreak18"D/L FlipBackArms,SwgFt,1/cs
BLS20FBD-ELR	BlueStreak 20"D/FlipBackArms,ElvLeg,1/ea
BLS20FBD-SF	BlueStreak 20"D/FlipBackArms,SwgFt,1/ea

SPECIFICATIONS

Description

Back Height	16"
Carton Shipping Weight	47 lbs
Height	36"
Length with Riggings	42"
Product Weight	41.25 lbs
Seat Dimensions with 18 Seat	18" (W) x 16" (D)
Seat Dimensions with 20 Seat	20" (W) x 16" (D)
Seat To Floor Height	19.5"
Warranty	Limited 3 Year
Weight Capacity	250 lbs
Width Closed	12.5"
Width Open	24" 26"

ACCESSORIES

Description

STDS833N	Rear Anti Tipper w.wheels SSP&CS,1/pr, PR
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA

ACCESSORIES	Description
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS831	Heel Loops, 1pr, BX
PARTS	Description
STDS4S2209	Back Upholstery 16",Blue Streak,4S,1ea, EA
STDS4S2208	Back Upholstery,BS, SSP1, 4S, 1/ea, EA
STDS4S2210	Back Upholstery 20",Blue Streak,4S,1ea, EA
STDS5N2202	Saddle Assy BlueStreak,2/pack, PKG
STDS3J2416	SeatUphol.16"Cruiser3,1/ea, EA
STDS3J2418	Seat Uphol Nylon 18, SSP1,Cruiser,1/ea, EA
STDS3J2420	SeatUphol.20"Cruiser3,1/ea, EA
STDS4S22X8	Cross Frame,Blue Streak 18", 4S,1/ea, EA
STDS4S2207	Side Frame,Left,Blue Streak,4S,1ea, EA
STDS5N2231	Caster Axle & Nut,BlueStreak,2 sets/pack, PKG
STDS4S22SF	Foot Rest Set,BlueStreak,4S,1set/pkg, PKG
STDS5N2201	Calf Strap BlueStreak, 1 ea/pkg, PKG
STDS5N2210	SlideTube&Bushing,BlueStreak,1/ea, EA
STDS5N22DL	L. Arm Assy. Complete BlueStreak,1/ea, EA
STDS5N22DR	R. Arm Assy. Complete BlueStreak,1/ea, EA
STDS4S2211	Left Side Guard W/Hrdwr,Blue Streak,1/ea, EA
STDS4S2212	Right Side Guard W/Hrdwr,Blue Streak,1ea, EA
STDS8M4603	BrakeTipSleeve,Cruiser,SSP, 1/ea, EA
STDS5N2219	Removable Arm Hardware, Blue Streak, 1ea, PKG
STDSDSKVB	UniversalWheelChairVinylDeskArmPad,1/ea, EA